2 × Taq Green Master Mix

Size:

Cat.#	Size
BG10109S	40 R(50 µL vol)
BG10109M	200 R(50 µL vol)
BG10109L	1000 R(50 µL vol)

Source: *E.coli* strain with a cloned *pol* gene from *Thermus aquaticus* YT1.

Storage temperature: -20°C

Exp.: 12 months(-15°C ~ -25°C)

Applications:

1. Routine PCR amplification of DNA fragments up to 5 kb.
2. PCR product for TA cloning.
3. Colony PCR.

website: syngentech.com

Description:

2 × Taq Green Master Mix contains Taq DNA Polymerase, dNTPs, other buffer components and stabilizers. Taq DNA Polymerase is a highly thermostable DNA polymerase, that possesses 5' → 3' polymerase activity, 5' → 3' exonuclease activity, but no detectable 3' → 5' exonuclease (proofreading) activity. Taq DNA Polymerase is extremely sensitive and able to amplify trace amounts of DNA in the sample reaction.

2 × Taq Green Master Mix is recommended for any amplification reaction that will be visualized by agarose gel electrophoresis and ethidium bromide staining. It is not recommended when doing any downstream applications in absorbance or fluorescence excitation, as the yellow and blue dyes in the reaction buffer may interfere with these applications.

Reagent supplied:

Cat.#	Size	Quantity	Nuclease-Free Water
BG10109S	40 R(50 µL vol)	1 mL	1 mL
BG10109M	200 R(50 µL vol)	5 × 1 mL	5 × 1 mL
BG10109L	1000 R(50 µL vol)	25 × 1 mL	25 × 1 mL

Protocol:

To prepare several parallel reactions with minimized pipetting errors, prepare a PCR master mix by mixing Nuclease-Free water, primers and 2 × Taq Master Mix. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube on ice and add the following component for each 25 µL or 50 µL reaction:

Component	25 µL reaction	50 µL reaction	Final concentration
10 µM forward primer	0.5 µL	1 µL	200 nM
10 µM reverse primer	0.5 µL	1 µL	200 nM
Template DNA	variable	variable	0.1-100 ng
2 × Taq Green Master Mix	12.5 µL	25 µL	1 ×
Nuclease-Free water	to 25 µL	to 50 µL	-

3. Gently vortex the samples and spin down.
4. Perform PCR using recommended thermal cycling conditions:

Step	Temperature	Time	Number of cycles
Initial Denaturation	94°C	2-5 min	1
Denaturation	94°C	15-30 s	25-40
Annealing	Tm-5°C	10-30 s	
Extension	68°C	10-60 s/kb	
Final Extension	68°C	5-15 min	1

General guidelines:

1. Template:

High quality, purified DNA templates greatly enhances the success of PCR reaction. Recommended amounts of DNA template for a 50 µL reaction are as follows:

DNA	Amount
Genomic DNA	1 ng-1 µg
Plasmid or Viral DNA	1 pg-1 µg

2. Primers:

Oligonucleotide primers are generally 20-40 nucleotides in length and ideally have a GC content of 40-60%. Software such as Oligo 7 can be used to design or analyze primers. The final concentration of each primer in a PCR reaction may be 0.05-1 µM, typically 0.1-0.5 µM.

3. Denaturation:

An initial denaturation of 2 min at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences or crude samples, a longer denaturation of 2-8 min at 94°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 2-5 min denaturation at 94°C is recommended.

4. Annealing:

The annealing step is typically 15-60 s. Annealing temperature is based on the Tm of the primer pair and is typically 45-68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting at 5°C below the calculated Tm.

5. Extension:

The recommended extension temperature is 68°C. Extension time are generally 20-30 s/kb for complex, genomic samples, but can be reduced to 10 s/kb for simple templates (plasmid,

virus). To determine specific extension time, we recommend 1 min/kb and reduce the extension time until 10 s/kb. 2 min/kb is recommended for crude samples. A final extension of 5 min at 68°C is recommended.

6. Cycle number:

Generally, 25-35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number templates.

7. PCR product:

The PCR products generated by Taq DNA Polymerase contain dA overhangs at the 3' end, therefore can be ligated to dT/dU-overhang vectors.

References:

1. Chien A, Edgar D B, Trela J M. Deoxyribonucleic-acid polymerase from extreme thermophile *Thermus-aquaticus*[J]. *Journal of Bacteriology*, 1976, 127.
2. Kaledin A S, Sliusarenko A G, Gorodetski S I. Isolation and properties of DNA polymerase from extreme thermophilic bacteria *Thermus aquaticus* YT-1[J]. *Biokhimiia*, 1980, 45(4): 644-651.

Research only, not for therapeutic or diagnostic purposes.

V1.0.3 © 13-14 Beijing Syngentech Co., Ltd. All rights reserved

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V1.0.3 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved