## Gene Taq™ DNA Polymerase with Green Reaction Buffer

### Size:

Cat.#	Size
BG10101S	100 units
BG10101M	500 units
BG10101L	2500 units

**Concentration:** 5 units/μL

**Source:** *E.coli* strain with a cloned pol gene from *Thermus aquaticus* YT1.

**Storage temperature:** -20°C

**Exp.:** 12 months(-15°C ~ -25°C)

### Applications:

1. Routine PCR amplification of DNA fragments up to 5 kb.
2. PCR product for TA cloning.
3. Colony PCR.

www.syngene.tech

### Description:

Gene Taq™ DNA polymerase is a highly thermostable DNA polymerase, that possesses 5' → 3' polymerase activity, 5' → 3' exonuclease activity, but no detectable 3' → 5' exonuclease (proofreading) activity. Gene Taq™ DNA Polymerase is extremely sensitive and able to amplify trace amounts of DNA in the sample reaction.

Gene Taq™ DNA polymerase with Green Reaction Buffer is recommended for any amplification reaction that will be visualized by agarose gel electrophoresis and ethidium bromide staining. Gene Taq™ DNA polymerase with Green Reaction Buffer is not recommended when doing any downstream applications in absorbance or fluorescence excitation, as the yellow and blue dyes in the reaction buffer may interfere with these applications.

### Reagent supplied:

Cat.#	Size	Reagent	Quantity
BG10101S	100 units	5xGene Taq™ Green Reaction Buffer	1.5 mL
		25 mM MgCl <sub>2</sub>	1.5 mL
BG10101M	500 units	5xGene Taq™ Green Reaction Buffer	4 x 1.5 mL
		25 mM MgCl <sub>2</sub>	3 x 1.5 mL
BG10101L	2500 units	5xGene Taq™ Green Reaction Buffer	20 x 1.5 mL
		25 mM MgCl <sub>2</sub>	15 x 1.5 mL

### Quality control:

1. Amplification Ability: 25 cycles of PCR amplification of 5 ng Lambda DNA in the expected 5 kb product.
2. Nuclease Activity: No contaminating endonuclease or exonuclease activity has been detected.
3. Protein Purity (SDS-PAGE): Gene Taq™ DNA Polymerase

has >99% pure as judged by SDS-PAGE.

4. qPCR DNA Contamination(*E.coli* Genomic): A minimum of 5μL of Gene Taq™ DNA Polymerase is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA. The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

### Definition of activity unit:

One unit of the enzyme is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 min at 75 °C .

### Storage buffer:

The enzyme is supplied in: 10 mM Tris-HCl (pH7.4@ 25°C ), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% (v/v) Tween 20, 0.5% (v/v) NP-40 and 50% (v/v) glycerol.

### Protocol:

To prepare several parallel reactions with minimized pipetting errors, prepare a PCR master mix by mixing Nuclease-Free water, buffer, Mg<sup>2+</sup>, dNTPs, primers and Gene Taq™ DNA polymerase. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube on ice and add the following component for each 25 μL or 50 μL reaction:

Component	25 μL reaction	50 μL reaction	Final concentration
5 × Gene Taq™ Reaction Buffer	5 μL	10 μL	1 ×
MgCl <sub>2</sub> solution, 25 mM	1-4 μL	2-8 μL	1-4 mM
dNTP Mix (10 mM each)	0.5 μL	1 μL	200 μM each
10 μM forward primer	0.5 μL	1 μL	200 nM
10 μM reverse primer	0.5 μL	1 μL	200 nM
Template DNA	variable	variable	0.1-100 ng
GeneTaq DNA Polymerase	0.125 μL	0.25 μL	–
Nuclease-Free water	to 25 μL	to 50 μL	–

3. Gently vortex the samples and spin down.

4. Perform PCR using recommended thermal cycling conditions:

Step	Temperature	Time	Number of cycles
Initial Denaturation	94°C	2-5 min	1
Denaturation	94°C	15-30 s	25-40
Annealing	Tm-5°C	10-30 s	
Extension	68°C	10-60 s/kb	
Final Extension	68°C	5-15 min	1

### General guidelines:

#### 1. Template

High quality, purified DNA templates greatly enhances the success of PCR reaction. Recommended amounts of DNA template for a 50 μL reaction are as follows:

DNA	Amount
Genomic DNA	1 ng-1 μg
Plasmid or Viral DNA	1 pg-1 μg

#### 2. Primers:

Oligonucleotide primers are generally 20-40 nucleotides in length and ideally have a GC content of 40-60%. Software such as Oligo 7 can be used to design or analyze primers. The final concentration of each primer in a PCR reaction may be 0.05-1 μM, typically 0.1-0.5 μM.

#### 3. Mg<sup>2+</sup> and additives:

Mg<sup>2+</sup> concentration of 1.5-2.0 mM is optimal for most PCR products generated with Gene Taq™ DNA Polymerase. The final Mg<sup>2+</sup> concentration in 1 × Gene Taq™ Reaction Buffer is 1.5 mM. This supports satisfactory amplification of most amplicons. However, Mg<sup>2+</sup> can be further optimized in 0.5 or 1.0 mM increments using MgCl<sub>2</sub>. Amplification of some difficult templates, like GC-rich sequences, may be improved with additives, such as DMSO or formamide.

#### 4. Denaturation:

An initial denaturation of 2 min at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences or crude samples, a longer denaturation of 2-8 min at 94°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 2-5 min denaturation at 94°C is recommended.

#### 5. Annealing:

The annealing step is typically 15-60 s. Annealing temperature is based on the Tm of the primer pair and is typically 45-68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting at 5°C below the calculated Tm.

## 6. Extension:

The recommended extension temperature is 68°C. Extension time are generally 20-30 s/kb for complex, genomic samples, but can be reduced to 10 s/kb for simple templates (plasmid, virus). To determine specific extension time, we recommend 1 min/kb and reduce the extension time until 10 s/kb. 2 min/kb is recommended for crude samples. A final extension of 5 min at 68°C is recommended.

## 7. Cycle number:

Generally, 25-35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number templates.

## 8. PCR product:

The PCR products generated by GeneTaq™ DNA Polymerase contain dA overhangs at the 3' end, therefore can be ligated to dT/dU-overhang vectors.

## References:

- Chien A, Edgar D B, Trela J M. Deoxyribonucleic-acid polymerase from extreme thermophile Thermus-aquaticus[J]. *Journal of Bacteriology*, 1976, 127.
- Kaledin A S, Sliusarenko A G, Gorodetski S I. Isolation and properties of DNA polymerase from extreme thermophilic bacteria Thermus aquaticus YT-1[J]. *Biokhimiia*, 1980, 45(4): 644-651.

Research only, not for therapeutic or diagnostic purposes.

V1.0.2 © 13-14 Beijing Syngentech Co., Ltd. All rights reserved

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