Quick Ligase™ Kit

Size:

Cat#.	Size
BG10202S	30 reactions
BG10202L	150 reactions

Storage temperature: -20°C .

Exp.: 12 months(-15°C ~ -25°C)

Applications:

1. Cloning of restriction fragments
2. Cloning of PCR products
3. Joining of double stranded linkers and adapters to DNA
4. Nick repair in duplex DNA, RNA or DNA/RNA hybrids

Description:

Quick Ligase™ catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. With the 2 × Quick ligase™ Buffer, the enzyme can ligate DNA sticky or blunt end fragments in 5 min at room temperature.

Reagent Supplied:

Cat#.	Size	Reagent	No.
BG10202S	30 reactions	2 × Quick Ligase™ Reaction Buffer	1 mL
BG10202L	150 reactions	2 × Quick Ligase™ Reaction Buffer	2 × 1 mL

Quality Control:

1. Protein Purity (SDS-PAGE): Quick Ligase™ is >99% pure as judged by SDS-PAGE.
2. Non-Specific DNase Activity (16 h): The product is tested for non-specific nuclease degradation in a reaction containing a DNA substrate. After incubation for 16 h there is no detectable degradation of the DNA substrate as determined by agarose gel electrophoresis.
3. Endonuclease Activity (Nicking): The product is tested in a reaction containing a supercoiled DNA substrate. After incubation for 4 h, < 5% DNA converted to the nicked form is determined by agarose gel electrophoresis.
4. qPCR for DNA Contamination(*E.coli* Genome): A minimum of 5 μL of Quick Ligase™ is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA. The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

Storage buffer:

10 mM Tris-HCl (pH7.5 @ 25°C), 50 mM KCl, 0.1 mM EDTA, 1 mM DTT and 50% (v/v) glycerol.

Protocol

1. Thaw the 2 × Quick DNA Ligase reaction buffer at room temperature and vortex vigorously to dissolve any precipitated material before use.
2. A molar ratio of 3:1 insert to vector(10-100ng) is recommended for the Reaction of DNA inserts to vectors. To calculate the amount of insert for a 3:1 ratio of vector to insert, use the formula:

$$X \text{ ng of insert} = \frac{3 \times \text{insert size (in bp)} \times \text{amount of vector (in ng)}}{\text{vector size (in bp)}}$$

Add the following reagents to an autoclaved, 0.2 mL microcentrifuge tube:

Reagent	Amount
2 × Quick Ligase™ Reaction Buffer	10 μL
Linearized vector DNA	10-100 ng
Insert DNA	variable
Quick Ligase™	1 μL
Nuclease-Free water	to 20 μL

3. Mix gently and centrifuge briefly.
4. Incubate at room temperature for 5 min.
5. Use 2-5 μL of the ligation mixture for transform or store at -20°C.
6. Do not heat inactivate. Heat inactivation dramatically reduces transformation efficiency.

Reference:

Remaut, E., Tsao, H., Fiers, W., Improved plasmid vectors with a thermoinducible expression and temperature-regulated runaway replication [J]. *Gene*, 1983, 22(1):103-113.

Research only, not for therapeutic or diagnostic purposes.

V1.0.2 © 13-14 Beijing Syngentech Co., Ltd. All rights reserved

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V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved