T4 DNA Ligase

Size:

Cat#.	Size
BG10201S	20000 U
BG10201L	100000 U

Concentration: 400 units/ μ L

Source: *E.coli* cells with a cloned T4 DNA ligase gene from bacteriophage T4.

Storage temperature: -20°C.

Exp.: 12 months(-15°C ~ -25°C)

Applications:

1. Cloning of restriction fragments
2. Cloning of PCR products
3. Joining of double stranded linkers and adapters to DNA
4. Nick repair in duplex DNA, RNA or DNA/RNA hybrids

www.syngentech.com

Description:

T4 DNA ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. With the 10 \times T4 DNA Ligase Buffer, the enzyme can ligate DNA sticky or blunt end fragments.

Reagent Supplied:

Cat#.	Size	Reagent	No.
BG10201S	20000 U	10 \times T4 DNA Ligase Reaction Buffer	1 mL
BG10201L	100000 U	10 \times T4 DNA Ligase Reaction Buffer	2 \times 1 mL

Quality Control:

1. Protein Purity (SDS-PAGE): T4 DNA Ligase has >99% pure judged by SDS-PAGE.
2. Non-Specific DNase Activity (16 h): The product is tested for non-specific nuclease degradation in a reaction containing a DNA substrate. After incubation for 16 h there is no detectable degradation of the DNA substrate as determined by agarose gel electrophoresis.
3. Endonuclease Activity (Nicking): The product is tested in a reaction containing a supercoiled DNA substrate. After incubation for 4 h, < 5% DNA converted to the nicked form is determined by agarose gel electrophoresis.
4. qPCR for DNA Contamination(*E.coli* Genome): A minimum of 5 μ L of T4 DNA Ligase is screened for the presence of *E.coli* genomic DNA using qPCR with Taqman Probe specific for *E.coli* 16S rRNA. The measured level of *E.coli* genomic DNA contamination is less than 10 copies of *E.coli* genome.

Definition of activity unit:

One unit of the enzyme is defined as the amount of enzyme required to reach 50% ligation of 100ng HindIII fragments of λ DNA in a total reaction volume of 20 μ L in 30 min at 16°C in 1 \times T4 DNA Ligase Buffer.

Storage buffer:

The enzyme is supplied in: 10 mM Tris-HCl (pH7.5 @ 25°C), 50 mM KCl, 0.1 mM EDTA, 1 mM DTT and 50% (v/v) glycerol.

Protocol:

1. Thaw the 10 \times T4 DNA Ligase Reaction Buffer at room temperature and vortex vigorously to dissolve any precipitated material before use.
2. A molar ratio of 3:1 insert: vector(10-100ng) is recommended for the Reaction of DNA inserts to vectors. To calculate the amount of insert for a 3:1 ratio of vector to insert, use the formula:

$$X \text{ ng of insert} = \frac{3 \times \text{insert size (in bp)} \times \text{amount of vector (in ng)}}{\text{vector size (in bp)}}$$

Add the following reagents to an autoclaved, 0.2 mL microcentrifuge tube:

Reagent	Amount
10 \times T4 DNA Ligase Reaction Buffer	2 μ L
Linearized vector DNA	10-100 ng
Insert DNA	variable
T4 DNA Ligase	1 μ L
Nuclease-Free water	to 20 μ L

3. Mix gently and centrifuge briefly.

4. Incubate at room temperature for 10 min or 16°C overnight for sticky ends. Incubate at room temperature for 2 h or 16°C overnight for blunt ends or single base overhang.
5. Use 2-5 μ L of the ligation mixture for transform or Store at -20°C.

Reference:

Remaut, E., Tsao, H., Fiers, W., Improved plasmid vectors with a thermoinducible expression and temperature-regulated runaway replication [J]. *Gene*, 1983, 22(1) :103-113.

Research only, not for therapeutic or diagnostic purposes.

V1.0.2 © 13-14 Beijing Syngentech Co., Ltd. All rights reserved

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V1.0.2 © 13-14 Beijing SyngenTech Co., Ltd. All rights reserved